

RITM

SCENAR®

PERSONAL DEVICES

www.scenar.com.au



*Non-invasive, drug free
therapy for pain relief*

RITMSCENAR® Sport D, Sport & Home


Sport D, Sport and Home, are electro-therapy devices with built-in electrodes that deliver complex electrical impulses through the skin to give pain relief.

SCENAR® interaction

The SCENAR® controls and automatically changes the impulse waveforms and frequencies to prevent the body from adapting to the stimulation.


RITMSCENAR® usage

RITMSCENAR® personal range of devices are designed for home users to manage their pain and to speed up the recovery process from acute and chronic pain, sport injury and post operative pain. The personal devices are light, compact and can be used almost anywhere.



"With consistent SCENAR treatment and maintenance, I was able to avoid surgery and retirement and return to competitive surfing without pain."

Layne Beachley, Australia
7x ASP Women's World
Surfing Champion






"Since 2009 SCENAR therapy is part of my daily injury prevention and management."

Nathan Cayless, Australia
Former NRL Parramatta
Eels Captain



"I am using RITM SCENAR for both injury and well-being while competing and travelling."

Torah Bright, Australia
World Snowboard Champion
and Olympic Gold medallist

FUNCTION/DEVICE	RITMSCENAR® Sport D	RITMSCENAR® Sport	RITMSCENAR® Home
			
Indication	LCD display	LED	LED
Skin contact indication	Yes	No	No
Dose Mode	Yes	Yes	Yes
Amplitude modulation	“3 : 1”	“3 : 1”	“3 : 1”
Frequency modulation	30Hz – 120Hz	30Hz – 120Hz	30Hz – 120Hz
Influencing frequencies	14Hz, 60Hz, 90Hz, 340Hz	14Hz, 60Hz, 90Hz, 340Hz	90Hz
Power supply	9V battery	9V battery	9V battery
Weight	0.2kg	0.2kg	0.2kg
Dimensions	13 x 5 x 2.5cm	13 x 5 x 2.5cm	13 x 5 x 2.5cm
Warranty	2 years	2 years	2 years

BRANCH OFFICE FOR ASIA PACIFIC REGION
RITM AUSTRALIA PTY LTD

7 Belsham Road Kariong NSW 2250 Australia
Tel: +61 2 8011 4217 • Fax: +61 2 9423 4681
info@scenar.com.au • www.scenar.com.au

DISTRIBUTOR



TRAINING

All RITMSCENAR® personal devices come with a comprehensive User's Manual, Flash training programme on USB Drive and a DVD that contains 19 practical videos showing device controls and the main treatment techniques.

Training seminars for personal device owners are also available on regular basis through out Australia.

For more information, please contact the Authorized Training Centre:

RITM SCENAR Institute Pty Ltd.
info@scenarinstitute.org
www.scenarinstitute.org



CERTIFICATES



Australia
ARTG
I40659

U.S.A.
FDA
K092117

Health
Canada
MDEL 4549