

DERMASWEEP

If we look at microdermabrasion systems versus the DermaSweep exfoliation, there are significant differences. If you are currently using a crystal system, you are dealing with the mess of crystals on the patient's skin as well as biohazard disposal of the crystals. Also crystal treatments utilize crystals that are all one size (usually about 100microns) which means that everyone receives the same level of treatment with only the vacuum pressure changing the treatment to make it less or more aggressive.

The different diamond head systems also offer multiple abrasion level without crystals, but the skin cells build up across the tip of the wands, which means you can lose efficacy as you progress through the treatment.

Microcurrent utilizes a low level of electricity. The electrical impulses trigger chemical reactions that are believed to affect collagen and elastin, provide circulatory benefits, and assist in product penetration through a process known as iontophoresis. Skin care products typically need to be Water-soluble and have a relatively small molecular size.

The DermaSweep offers multiple exfoliating bristles, which will give you the greatest treatment range for different skin types. These bristle tips offer you and your patient a very clean procedure. It is also very visual for the patient, as the filter disc will allow them to see the skin cells that are exfoliated. The bristle tips are exclusively designed to exfoliate while allowing for airflow through the bristles, which means you are able to take full advantage of the therapeutic power of the vacuum. Most systems make it difficult to properly use the vacuum to the fullest efficacy. The bristles truly make this a very different procedure which, from a hands-on perspective, really puts it in a class of its' own. The skin specific serums are then driven to the dermis through the Epi-Infusion process, which utilizes a silk bristle tip that drives solutions deep. Each skin specific solution is infused through an IV bag system and includes options for acne, dehydration, antioxidant (vitamin C), skin freshening and anti-aging.

The most remarkable difference is in patient reaction and efficacy. Patients are getting noticeable improvement with one treatment. Obviously like other treatments the ultimate benefits are seen after a series of treatments.

At first glance, it might appear as microderm and microcurrent offer the same type of treatment as the DermaSweep, the concept maybe the same, but the patient outcome and efficacy is completely different. The skin care practitioner is the primary factor in product penetration based strength of treatment and number of Epi-Infused passes you make over the skin. The DermaSweep involves a physical system for penetration, where you are in control of the process.