

GO-GO
ELITE TRAVELLER®

THE QUICKEST WAY
TO TRAVEL.

GO-GO ELITE TRAVELLER®

The Go-Go Elite Traveller®'s compact design allows it to easily manoeuvre in tight spaces while providing stable outdoor performance. Take the guesswork out of travel with simple, one hand Feather-touch disassembly. Go where you want to go, easily, with the Go-Go Elite Traveller.



Pride
Mobility Products
Australia Pty. Ltd.

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

Go-Go Elite Traveller®	3-wheel	4-wheel
Weight capacity:	125 kg (275 lbs.)	125 kg (275 lbs.)
Maximum speed ² :	Up to 7 km/h (4.25 mph)	Up to 7 km/h (4.25 mph)
Ground clearance ¹ :	6.5 cm (2.625") at midframe	6.35 cm (2.5") at midframe
Turning radius ¹ :	82.5 cm (32.5")	110 cm (43.25")
Overall length ¹ :	92.75 cm (36.5")	101 cm (39.75")
Overall width ¹ :	49 cm (19.25")	48.5 cm (19.125")
Tyres:		
Front:	5 cm x 20 cm (2" x 8") solid	5 cm x 18 cm (2" x 7") solid
Rear:	6.35 cm x 20 cm (2.5" x 8") solid	6.35 cm x 20 cm (2.5" x 8") solid
Anti-tips:	Rear	Rear
Transport weight:		
Total weight w/o battery pack:	34 kg (75 lbs.)	35.5 kg (78.5 lbs.)
Heaviest piece when disassembled:	12.75 kg (28 lbs.) (front section)	15 kg (33.5 lbs.) (front section)
Drive System:	24-volt DC motor, sealed mini-transaxle	
Dual Braking System:	Electronic, regenerative, and electromechanical	
Standard seating:	Foldable, compact, vinyl-covered molded plastic	
Battery Charger:	Off-board, 2-amp	
Per-charge range ² :	Up to 22.5 km (14 miles) (3 or 4-wheel w/ 2 x 18 AH) standard	
Battery Requirements ³ :	Two 12-volt, deep cycle (standard) Size: 2 x 18 AH Weight: 18 AH: 15 kg (33 lbs.) battery pack	

Standard Colours:



Blue



Red



Silver

3 sets of colours with every scooter!

Features

- Feather-touch disassembly permits simple frame separation with only one hand
- Easy drive tiller with wraparound handles
- Pride's exclusive black non-scuffing tyres with tri-spoke wheels
- Auto-connecting front to rear cable
- Includes 3 sets of easily changeable coloured panels in 3 exciting colours - Red, Blue and Silver
- Front frame-mounted seat post offers maximum stability
- All-new frame design easily disassembles into 5 super light-weight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket included
- Removable deck basket for additional storage (4-wheel only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board
- 2 x 18 AH battery pack for maximum range per charge (standard)

Options

- Cup holder
- 2 x 12 AH battery pack for maximum portability



Feather-touch disassembly

¹ Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

² Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

³ Battery weight may vary based on manufacturer.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

21 Healey Rd., Dandenong, Victoria 3175
(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

ACN 088 609 661 ABN 66 088 609 661



Pride
Mobility Products
Australia Pty. Ltd.

Pathrider 10 Deluxe

SUPREME PERFORMANCE,
EXHILERATING STYLE,
AND EXCELLENT VALUE.

PATHRIDER 10 DELUXE

The sleek, sporty Pathrider 10 Deluxe scooter delivers smooth performance and all-new features. With a full complement of advanced standard features including, rugged suspension, high-visibility automatic front and rear lighting system, Feather-touch disassembly, LED kerb light, and wraparound easy-drive tiller, the Pathrider 10 Deluxe is a greater value than ever before.



Pride
Mobility Products
Australia Pty. Ltd.

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

Pathrider 10 Deluxe

Weight capacity:	159 kg (350 lbs.)
Maximum speed ¹ :	Up to 8.5 km/h (5.28 mph)
Ground clearance ² :	3.5 cm (1.38")
Turning radius:	137.16 cm (54")
Overall length ² :	118.75 cm (46.75")
Overall width ² :	56.5 cm (22.24")
Front tyres:	9 cm x 25 cm (3.5" x 10")
Rear tyres:	10 cm x 25 cm (4" x 10")
Suspension:	Front and rear
Range per charge ¹ :	Up to 33.38 km (20.74 miles) with 40 AH batteries
Maximum obstacle climbing ability:	5 cm (2")
Total Weight Without Batteries:	84 kg (185 lbs.)
Heaviest piece when disassembled:	29.48 kg (65 lbs.)
Standard seating:	Type: CRS Junior seat Weight: 19 kg (42 lbs.) Material: Black Vinyl Dimensions: width – 44.5 cm (17.5") depth – 43 cm (17")
Drive System:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System:	Electronic, regenerative, and electromechanical
Battery Requirements ³ :	Two 12-volt, deep cycle Size: 35 AH or optional 40 AH Weight ⁴ : 35 AH: 11 kg (24.5 lbs.) each; 40 AH: 15 kg (33 lbs.) each
Battery Charger	5-amp, off-board

Colours:



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

⁴Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

Features

- Front and rear suspension
- High-visibility automatic front and rear lighting system
- Exclusive one-hand Feather-touch disassembly provides the ultimate quick transport system
- Long-lasting, high-intensity LED headlight, tail light and indicators
- Pride's exclusive stylish, light weight, non-marking, low profile wheels
- Auto-connecting front-to-rear harness
- Easy-drive tiller with wraparound handles
- New back-lit battery meter for increased accuracy and visibility in any setting
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger
- Rear view mirror standard

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Pathrider™ 10

STYLISH DESIGN & PREMIUM
PERFORMANCE AT AN
EXCELLENT VALUE.

PATHRIDER™ 10

The sleek, sporty **Pathrider™ 10** scooter delivers high performance operation, all-new features and feather-touch disassembly. With a full complement of advanced standard features exclusive low-profile, puncture-proof tyres, back-lit battery gauge, and wraparound easy drive tiller, the Pathrider 10 is a greater value than ever before.



Pride®
Mobility Products
Australia Pty. Ltd.


03 9706 4611 (Tel) ■ 03 9706 4622 (Fax) ■ www.pridemobility.com

Pathrider™ 10

4-wheel scooter

Weight capacity:	159 kg (350 lbs.)
Maximum speed ¹ :	Up to 8.5 km/h (5 mph)
Ground clearance ² :	5 cm (2")
Turning radius:	157.5 cm (62")
Overall length ² :	119 cm (47")
Overall width ² :	56.5 cm (22")
Tyre size:	Front & Rear: 9 cm x 25 cm (3.5" x 10") (puncture-proof)
Range per charge ¹ :	Up to 40 km (25 miles)
Maximum obstacle climbing ability:	5 cm (2")
Total Weight Without Batteries:	63.5 kg (140 lbs.)
Heaviest piece when disassembled:	32 kg (71 lbs.)
Standard seating:	Type: Foldable molded plastic Weight: 13 kg (28 lbs.) Material: Black Vinyl Dimensions: 46 cm (18") width 43 cm (17") depth
Drive System:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System:	Electronic, regenerative, and electromechanical
Battery Requirements ³ :	Two 12-volt, deep cycle Size: 35 AH (standard) or 40 AH (optional) Weight: 35 AH: 11 kg (24 lbs.) each; 40 AH: 15 kg (32.5 lbs.) each
Battery Charger:	3-amp, off-board

Colour:



MR MOBILITY P/L
THE SCOOTER MAN
163 SUNSHINE RD WEST FOOTSCRAY 3012
PH FAX 9314 9559 MOB 0407 056 960
www.mrmobility.com.au
ABN 39 107 967 724

Features

- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Easy drive tiller with wrap-around handles
- Pride's exclusive stylish, light weight, low profile, non-scuffing, puncture-proof wheels
- New back-lit LED battery meter for increased accuracy and visibility in any setting
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight, molded plastic seat with viscoelastic foam inserts for increased comfort
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or Gel-Cell type required.

NOTE: All specifications subject to change without notice.

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".

21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com
ACN 088 609 661 ABN 66 088 609 661

Tel: 03 9706 4611
Fax: 03 9706 4622



Pride[®]
Mobility Products
Australia Pty. Ltd.